

FGF11 RABBIT PAB

Cat.#: N225191

Product Name: Anti-FGF11 Rabbit pAb

Synonyms: FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3

UNIPROT ID: Q92914

Background: Probably involved in nervous system development and function.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human FGF11. AA range:21-70

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 25 kDa; Observed MW: 25 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
