

## FKBP51 RABBIT MAB

**Cat.#:** N262223

**Product Name:** Anti-FKBP51 Rabbit Monoclonal Antibody

**Synonyms:** P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10

**UNIPROT ID:** Q13451

**Background:** Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

**Immunogen:** A synthetic peptide of human FKBP51

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R08-8E5

**MW:** Calculated MW: 51 kDa; Observed MW: 51 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Hamster


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.