

FKBP51 RABBIT MAB

Cat.#: N262223

Product Name: Anti-FKBP51 Rabbit Monoclonal Antibody

Synonyms: P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10

UNIPROT ID: Q13451

Background: Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Immunogen: A synthetic peptide of human FKBP51

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-8E5

MW: Calculated MW: 51 kDa; Observed MW: 51 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Hamster


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.