

FLOTILLIN 1 RABBIT MAB

Cat.#: N262226

Product Name: Anti-Flotillin 1 Rabbit Monoclonal Antibody

Synonyms: Flotillin-1

UNIPROT ID: O75955

Background: May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Immunogen: A synthetic peptide of human Flotillin 1

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-5E9

MW: Calculated MW: 47 kDa; Observed MW: 47 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Flotillin 1 (green) in HeLa using Flotillin 1 antibody, and DAPI (blue)



Western blot analysis of Flotillin 1 in K562, HeLa lysates using Flotillin 1 antibody.