

GAD67 RABBIT MAB

Cat.#: N261670

Product Name: Anti-GAD67 Rabbit Monoclonal Antibody

Synonyms: Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67

UNIPROT ID: Q99259

Background: Catalyzes the production of GABA.

Immunogen: A synthetic peptide of human GAD67

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-2D9

MW: Calculated MW: 67 kDa; Observed MW: 67 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience GABAergic Neurons

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GAD67 in HeLa lysates using GAD67 antibody.