

## GAD67 RABBIT PAB

**Cat.#:** N225801

**Product Name:** Anti-GAD67 Rabbit pAb

**Synonyms:** Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67

**UNIPROT ID:** Q99259

**Background:** Catalyzes the production of GABA.

**Immunogen:** A synthesized peptide derived from human GAD67

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 67 kDa; Observed MW: 67 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience GABAergic Neurons

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GAD67 in (1) Jurkat lysates; (2) C6 lysates using GAD67 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---