

GDF6 RABBIT PAB

Cat.#: N225076

Product Name: Anti-GDF6 Rabbit pAb

Synonyms: GDF6; GDF16; Growth/differentiation factor 6; GDF-6; Growth/differentiation factor 16

UNIPROT ID: Q6KF10

Background: Growth factor that controls proliferation and cellular differentiation in the retina and bone formation.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human GDF6. AA range:311-360

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 51 kDa; Observed MW: 51 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GDF6 in 293T lysates using GDF6 antibody.



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.