

GEMIN7 RABBIT PAB

Cat.#: S219199

Product Name: Anti-GEMIN7 Rabbit Polyclonal Antibody

Synonyms: SIP3

UNIPROT ID: Q9H840 (Gene Accession - BC007793)

Background:

Immunogen: Fusion protein of human GEMIN7

Applications: ELISA, IHC

Recommended Dilutions: IHC: x which is required for pre-mRNA splicing in the nucleus. The encoded protein is found in the nucleoplasm;

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG

Purification: Antigen affinity purification


Species Reactivity: Human, Mouse

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing

 and in the cytoplasm. Three transcript variants encoding the same  protein have been found for this gene."
<http://store.sangon.com/productImage/PIC/219199/219199-IHC-1-P.png>

 The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using 219199(GEMIN7 Antibody) at dilution 1/50, on the right is treated with fusion protein.