

GRB2 RABBIT PAB

Cat.#: N225912

Product Name: Anti-GRB2 Rabbit pAb

Synonyms: GRB2; ASH; Growth factor receptor-bound protein 2; Adapter protein GRB2; Protein Ash; SH2/SH3 adapter GRB2

UNIPROT ID: P62993

Background: Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

Immunogen: A synthesized peptide derived from human GRB2

Applications: WB,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 25 kDa; Observed MW: 25 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GRB2 in HEK293 lysates using GRB2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
