

## GRK3 RABBIT PAB

**Cat.#:** N225956

**Product Name:** Anti-GRK3 Rabbit pAb

**Synonyms:** ADRBK2; BARK2; GRK3

**UNIPROT ID:** P35626

**Background:** Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

**Immunogen:** A synthesized peptide derived from human GRK3

**Applications:** WB,IHC-P,ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 80 kDa; Observed MW: 80 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GRK3 in Ramos lysates using GRK3 antibody.