

GRK3 RABBIT PAB

Cat.#: N226028

Product Name: Anti-GRK3 Rabbit pAb

Synonyms: ADRBK2; BARK2; GRK3

UNIPROT ID: P35626

Background: Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

Immunogen: A synthetic peptide of human GRK3

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 80 kDa; Observed MW: 80 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GRK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GRK3 in Jurkat lysates using GRK3 antibody.