

GST TAG MAB

Anti-GST Tag Mouse Monoclonal Antibody

Cat. #: 28005

Size: 100 µL

Symbol: GST-tag

Description: Anti-GST Tag Mouse Monoclonal Antibody

Background: The His-tag antibody recognizes any protein containing continuous 6 or more histidine residues.

Immunogen: GST

Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA 1:2000–1:5000

WB 1:1000–1:5000

IF 1:100–1:200

Concentration: 1 mg/ml

Host Species: Mouse

Format: Liquid

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites


Preservative: No

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes GST protein.

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



Western blot: Recombinant GST tag fusion protein was diluted with GST mouse monoclonal antibody (1;1000), and the protein loading amount was 50ng, and Western blot analysis was performed.