

GST3 RABBIT MAB

Cat.#: N262309

Product Name: Anti-GST3 Rabbit Monoclonal Antibody

Synonyms: GSTP1; FAES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1

UNIPROT ID: P09211

Background: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Immunogen: A synthetic peptide of human GST3/GST pi

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-2F4

MW: Calculated MW: 23 kDa; Observed MW: 23 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using GST3/GST pi antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GST3/GST pi in HepG2, HeLa lysates using GST3 antibody.