

HARS RABBIT MAB

Cat.#: N262313

Product Name: Anti-HARS Rabbit Monoclonal Antibody

Synonyms: HRS; CMT2W; USH3B

UNIPROT ID: P12081

Background: Cytoplasmic histidine--tRNA ligase (Probable). Plays a role in axon guidance.

Immunogen: A synthetic peptide of human HARS

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-1D2

MW: Calculated MW: 57 kDa; Observed MW: 57 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HARS in K562, HeLa lysates using HARS antibody.