

## HEME OXYGENASE 1 RABBIT MAB

**Cat.#:** N261674

**Product Name:** Anti-Heme Oxygenase 1 Rabbit Monoclonal Antibody

**Synonyms:** HO-1; HSP32; HMOX1; ; HO1

**UNIPROT ID:** P09601

**Background:** Hemeoxygenase (HO) is the rate-limiting enzyme in the catabolism of heme that results in the release of carbon monoxide, iron, and biliverdin. The products of this enzymatic reaction play important biological roles in antioxidant, anti-inflammatory and cytoprotective functions. Hemeoxygenase comprises two isozymes, including the constitutively expressed HO-2 isozyme and the inducible HO-1 isozyme.

**Immunogen:** A synthetic peptide of human Heme Oxygenase 1

**Applications:** WB,ICC/IF,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R03-8A2

**MW:** Calculated MW: 33 kDa; Observed MW: 33 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Mouse


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Heme Oxygenase 1 in 3T3 lysates using Heme Oxygenase 1 antibody.



Immunocytochemistry analysis of Heme Oxygenase 1 (green) in U87-MG using Heme Oxygenase 1 antibody, and DAPI (blue).