

HGS RABBIT PAB

Cat.#: N225173

Product Name: Anti-HGS Rabbit pAb

Synonyms: HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110

UNIPROT ID: O14964

Background: Involved in intracellular signal transduction mediated by cytokines and growth factors.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human HGS. AA range:511-560

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 86 kDa; Observed MW: 86 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HGS in mouse eye, mouse kidney, mouse testis lysates using HGS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
