

## HKDC1 RABBIT MAB

**Cat.#:** N263367

**Product Name:** Anti-HKDC1 Rabbit Monoclonal Antibody

**Synonyms:** Hexokinase domain-containing protein 1

**UNIPROT ID:** Q2TB90

**Background:** Catalyzes the phosphorylation of hexose to hexose 6-phosphate, although at very low level compared to other hexokinases.

**Immunogen:** Recombinant protein of human HKDC1

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R06-814

**MW:** Calculated MW: 103 kDa; Observed MW: 103 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HKDC1 in Hela, CHO-K1, Ramos lysates using HKDC1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---