

HNF 4 ALPHA RABBIT MAB

Cat.#: N263371

Product Name: Anti-HNF 4 alpha Rabbit Monoclonal Antibody

Synonyms: TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7; HNF4a8; HNF4a9; NR2A21; HNF4alpha

UNIPROT ID: P41235

Background: HNF-4-alpha is a transcription factor that belongs to the steroid hormone receptor superfamily. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha.

Immunogen: Recombinant protein of human HNF-4-alpha

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-4F0

MW: Calculated MW: 53 kDa; Observed MW: 53 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HNF4alpha in A549, MCF-7 lysates using HNF 4 alpha antibody.