

## HNF 4 ALPHA RABBIT MAB

**Cat.#:** N262345

**Product Name:** Anti-HNF 4 alpha Rabbit Monoclonal Antibody

**Synonyms:** Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha

**UNIPROT ID:** P49698

**Background:** Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine. Miscellaneous Binds fatty acids.

**Immunogen:** Recombinant protein of mouse HNF-4-alpha

**Applications:** WB, IHC-F, IHC-P, ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R02-7J6

**MW:** Calculated MW: 53 kDa; Observed MW: 53 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Mouse


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cardiovascular

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HNF4alpha in Raw264.7 lysates using HNF 4 alpha antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---