

HRPT2 RABBIT PAB

Cat.#: N225226

Product Name: Anti-HRPT2 Rabbit pAb

Synonyms: Parafibromin (Cell division cycle protein 73 homolog; Hyperparathyroidism 2 protein)

UNIPROT ID: Q6P1J9

Background: Tumor suppressor probably involved in transcriptional and post-transcriptional control pathways.

Immunogen: Synthesized peptide derived from human Parafibromin. at AA range: 51-100

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 61 kDa; Observed MW: 61 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HRPT2 in mouse liver, mouse brain lysates using HRPT2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
