

## HRPT2 RABBIT PAB

**Cat.#:** N225226

**Product Name:** Anti-HRPT2 Rabbit pAb

**Synonyms:** Parafibromin (Cell division cycle protein 73 homolog; Hyperparathyroidism 2 protein)

**UNIPROT ID:** Q6P1J9

**Background:** Tumor suppressor probably involved in transcriptional and post-transcriptional control pathways.

**Immunogen:** Synthesized peptide derived from human Parafibromin. at AA range: 51-100

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 61 kDa; Observed MW: 61 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HRPT2 in mouse liver, mouse brain lysates using HRPT2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---