

LGALS14 RABBIT MAB

Cat.#: N262468

Product Name: Anti-LGALS14 Rabbit Monoclonal Antibody

Synonyms: Charcot-Leyden crystal protein 2; CLC2; Galectin-14; Gal-14

UNIPROT ID: Q8TCE9

Background: Binds beta-galactoside and lactose.

Immunogen: A synthetic peptide of human LGALS14

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-6C1

MW: Calculated MW: 16 kDa; Observed MW: 16 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LGALS14 in C6 lysates using LGALS14 antibody.



Immunocytochemistry analysis of LGALS14 (green) in K562 using LGALS14 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
