

## LILRB1 RABBIT MAB

**Cat.#:** N263404

**Product Name:** Anti-LILRB1 Rabbit Monoclonal Antibody

**Synonyms:** ILT2; LIR1; MIR7; PIRB; CD85J; ILT-2; LIR-1; MIR-7; PIR-B

**UNIPROT ID:** Q8NHL6

**Background:** Receptor for class I MHC antigens.

**Immunogen:** Recombinant protein of human LILRB1

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R05-2A4

**MW:** Calculated MW: 71 kDa; Observed MW: 110 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.