

LILRB1 RABBIT MAB

Cat.#: N263404

Product Name: Anti-LILRB1 Rabbit Monoclonal Antibody

Synonyms: ILT2; LIR1; MIR7; PIRB; CD85J; ILT-2; LIR-1; MIR-7; PIR-B

UNIPROT ID: Q8NHL6

Background: Receptor for class I MHC antigens.

Immunogen: Recombinant protein of human LILRB1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-2A4

MW: Calculated MW: 71 kDa; Observed MW: 110 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.