

## LLGL1 RABBIT MAB

**Cat.#:** N263407

**Product Name:** Anti-LLGL1 Rabbit Monoclonal Antibody

**Synonyms:** DLG4; HUGL; LLGL; Lgll; MglI; HUGL1; HUGL-1

**UNIPROT ID:** Q15334

**Background:** Cortical cytoskeleton protein found in a complex involved in maintaining cell polarity and epithelial integrity.

**Immunogen:** Recombinant protein of human LLGL1

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R05-0B8

**MW:** Calculated MW: 115 kDa; Observed MW: 115 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LLGL1 in A549, U2OS lysates using LLGL1 antibody.