

LRG1 RABBIT MAB

Cat.#: N262477

Product Name: Anti-LRG1 Rabbit Monoclonal Antibody

Synonyms: LRG1; LRG; Leucine-rich alpha-2-glycoprotein; LRG

UNIPROT ID: P02750

Background: Transforming growth factor beta receptor binding.

Immunogen: A synthetic peptide of human LRG1

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-413

MW: Calculated MW: 38 kDa; Observed MW: 48 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LRG1 in Human serum lysates using LRG1 antibody.