

MKI67 RABBIT PAB

Cat.#: S219509

Product Name: Anti-MKI67 Rabbit Polyclonal Antibody

Synonyms: KIA; MIB-; MIB-1; PPP1R105

UNIPROT ID: P46013 (Gene Accession - NP_002408)

Background: This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

Immunogen: Fusion protein of human MKI67

Applications: ELISA, IHC

Recommended Dilutions: IHC: 100-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cancer, Neuroscience


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human tonsil tissue using 219509(MKI67 Antibody) at a dilution of 1/85(Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with the fusion protein and then with 219509(Anti-MKI67 Antibody) at dilution 1/85.



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using 219509(Anti-MKI67 Antibody) at a dilution of 1/85.



In comparison with the IHC on the left, the same paraffin-embedded Human breast cancer tissue is first treated with fusion protein and then with D227440(Anti-MKI67 Antibody) at dilution 1/85.