

MKLP1 RABBIT PAB

Cat.#: N225552

Product Name: Anti-MKLP1 Rabbit pAb

Synonyms: CHO1; KNSL5; MKLP1; MKLP-1

UNIPROT ID: Q02241

Background: Plus-end-directed motor enzyme that moves antiparallel microtubules in vitro. Localizes to the interzone of mitotic spindles.

Immunogen: A synthesized peptide derived from human MKLP1

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 110 kDa; Observed MW: 98 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human, Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MKLP1 in HeLa lysates using MKLP1 antibody.