


**MMAE(11B2) RABBIT MAB****Cat.#:** 28573**Product Name:** Anti-MMAE(11B2) Rabbit Monoclonal Antibody**Synonyms:** N.A.**Description:** Anti-MMAE antibody(11B2) Rabbit Monoclonal Antibody**Background:** N.A.**Applications:** ELISA**Recommended Dilutions:** ELISA 1/5000-10000**Host Species:** Rabbit**Isotype:** Rabbit IgG**Purification:** Purified from cell culture supernatant by affinity chromatography**Species Reactivity:** N.A. MMAE**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2µg/per well). Serial diluted anti-MMAE monoclonal antibody (28573) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-MMAE monoclonal antibody binding with IgG-MMAE is 2.538ng/ml.

