

## MMAE(8B4) RABBIT MAB

**Cat.#:** 28570

**Product Name:** Anti-MMAE(8B4) Rabbit Monoclonal Antibody

**Synonyms:** N.A.

**Description:** Anti-MMAE antibody(8B4) Rabbit Monoclonal Antibody

**Background:** N.A.

**Applications:** ELISA

**Recommended Dilutions:** ELISA 1/5000-10000

**Host Species:** Rabbit

**Isotype:** Rabbit IgG

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** N.A. MMAE

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).

ELISA assay to evaluate Anti-MMAE antibody  
0.2 $\mu\text{g}$  Human IgG-MMAE per well




Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2 $\mu\text{g}$ /per well). Serial diluted anti-MMAE monoclonal antibody (28570) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC<sub>50</sub> for anti-MMAE monoclonal antibody binding with IgG-MMAE is 2.274ng/ml.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---