

MMP9 RABBIT MAB

Cat.#: N262530

Product Name: Anti-MMP9 Rabbit Monoclonal Antibody

Synonyms: 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB

UNIPROT ID: P41245

Background: Matrix metalloproteinase that plays an essential role in local proteolysis of the extracellular matrix and in leukocyte migration (By similarity).

Immunogen: Recombinant protein of mouse MMP9

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-3A5

MW: Calculated MW: 81 kDa; Observed MW: 92 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MMP9 in rat lung lysates using MMP9 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
