

MTTFA RABBIT PAB

Cat.#: N225918

Product Name: Anti-mtTFA Rabbit pAb

Synonyms: MtTF1; mtTFA; TCF6; TCF6L1; TCF6L2; TCF6L3; TFAM; Transcription factor 6

UNIPROT ID: Q00059

Background: Binds to the mitochondrial light strand promoter and functions in mitochondrial transcription regulation. Required for accurate and efficient promoter recognition by the mitochondrial RNA polymerase.

Immunogen: A synthesized peptide derived from human mtTFA

Applications: WB,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 29 kDa; Observed MW: 24 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of mtTFA in MCF-7 lysates using mtTFA antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
