

MYC TAG MAB

Anti-MYC Tag Mouse Monoclonal Antibody

Cat. #: 28001

Size: 100 μ L

Symbol: C-Myc

Description: Anti-MYC Tag Mouse Monoclonal Antibody

Background: The Anti-MYC antibody recognizes any protein containing Myc Tag.

Immunogen: CEQKLISEEDL

Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA 1:2000-1:5000

WB 1:1000-1:5000

IF 1:100-1:200

Concentration: 1 mg/ml

Host Species: Mouse

Format: Liquid

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Preservative: No


Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes protein contain Myc .

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

c-Mycprotein



Western blot: Recombinant c-Myc tag fusion protein was diluted with c-Myc mouse monoclonal antibody (1:1000), and the protein loading amount was 50 ng, and Western blot analysis was performed.