

MYC TAG MAB**Anti-MYC Tag Mouse Monoclonal Antibody****Cat. #:** 28001**Size:** 100 µL**Symbol:** C-Myc**Description:** Anti-MYC Tag Mouse Monoclonal Antibody**Background:** The Anti-MYC antibody recognizes any protein containing Myc Tag.**Immunogen:** CEQKLISEEDL**Applications:** ELISA, WB, IF, IHC**Recommended Dilutions:**

ELISA 1:2000–1:5000


WB 1:1000–1:5000

IF 1:100–1:200

Concentration: 1 mg/ml**Host Species:** Mouse**Format:** Liquid**Clonality:** Monoclonal**Isotype:** IgG**Purity:** Purified from ascites**Preservative:** No**Constituents:** PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol**Species Reactivity:** Recognizes protein contain Myc .**Storage Conditions:** Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

c-Mycprotein



Western blot: Recombinant c-Myc tag fusion protein was diluted with c-Myc mouse monoclonal antibody (1:1000), and the protein loading amount was 50 ng, and Western blot analysis was performed.