

## N4BP1 RABBIT MAB

**Cat.#:** N262567

**Product Name:** Anti-N4BP1 Rabbit Monoclonal Antibody

**Synonyms:** C81621; A1481586

**UNIPROT ID:** Q6A037

**Background:** Potent suppressor of cytokine production that acts as a regulator of innate immune signaling and inflammation.

**Immunogen:** Recombinant protein of mouse N4BP1

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R08-6E2

**MW:** Calculated MW: 99 kDa; Observed MW: 99 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of N4BP1 in K562, HeLa lysates using N4BP1 antibody.