

NECTIN 2 RABBIT MAB

Cat.#: N263446

Product Name: Anti-Nectin 2 Rabbit Monoclonal Antibody

Synonyms: HVEB; PRR2; CD112; PVRR2

UNIPROT ID: Q92692

Background: Probable cell adhesion protein.

Immunogen: A synthetic peptide of human Nectin 2

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-8F3

MW: Calculated MW: 58 kDa; Observed MW: 70-80 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Nectin 2 (green) in HeLa using Nectin 2 antibody, and DAPI (blue)



Western blot analysis of Nectin 2 in HeLa, A549, HL-60, U251, U87-MG lysates using Nectin 2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
