

NEDD4 RABBIT MAB

Cat.#: N263125

Product Name: Anti-NEDD4 Rabbit Monoclonal Antibody

Synonyms: RPF1; NEDD4-1

UNIPROT ID: P46934

Background: E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Involved in the pathway leading to the degradation of VEGFR-2/KDFR, independently of its ubiquitin-ligase activity.

Immunogen: A synthetic peptide of human NEDD4

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-7C3

MW: Calculated MW: 149 kDa; Observed MW: 115 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NEDD4 in HeLa, A549, HL-60, U2OS, U87-MG lysates using NEDD4 antibody.



Western blot analysis of NEDD4 in K562, HCT116, 3T3, C6 lysates using NEDD4 antibody.