

## NEUROMEDIN B RECEPTOR RABBIT PAB

**Cat.#:** N225532

**Product Name:** Anti-Neuromedin B Receptor Rabbit pAb

**Synonyms:** Neuromedin B receptor antibody; Bombesin Receptor 1 antibody; BBI antibody

**UNIPROT ID:** P28336

**Background:** Receptor for neuromedin-B.

**Immunogen:** Synthetic Peptide of NMBR

**Applications:** WB,IHC-P

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 43 kDa; Observed MW: 50 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience yyyyyy Neuromedin B Receptor

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Neuromedin B Receptor in Human Brain, rat Brain lysates using NMBR(Bombesin Receptor) antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---