

## NOXIN RABBIT PAB

**Cat.#:** N225194

**Product Name:** Anti-NOXIN Rabbit pAb

**Synonyms:** NOXIN; C11orf82; Nitric oxide-inducible gene protein

**UNIPROT ID:** Q8IXT1

**Background:** May be an anti-apoptotic protein involved in DNA repair or cell survival.

**Immunogen:** The antiserum was produced against synthesized peptide derived from the Internal region of human NOXIN. AA range:111-160

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 112 kDa; Observed MW: 85 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NOXIN in 823, 293T, AD293, HeLa, HepG2 lysates using NOXIN antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---