

NUDT9 RABBIT MAB

Cat.#: N262633

Product Name: Anti-NUDT9 Rabbit Monoclonal Antibody

Synonyms: NUDT10

UNIPROT ID: Q9BW91

Background: Hydrolyzes ADP-ribose (ADPR) to AMP and ribose 5'-phosphate.

Immunogen: A synthetic peptide of human NUDT9

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-9C3

MW: Calculated MW: 39 kDa; Observed MW: 39 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of NUDT9 (green) in HeLa using NUDT9 antibody, and DAPI (blue)



Western blot analysis of NUDT9 in K562, HeLa lysates using NUDT9 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
