

NUP50 RABBIT MAB

Cat.#: N263465

Product Name: Anti-NUP50 Rabbit Monoclonal Antibody

Synonyms: NPAP60; NPAP60L

UNIPROT ID: Q9UKX7

Background: Component of the nuclear pore complex that has a direct role in nuclear protein import.

Immunogen: A synthetic peptide of human NUP50

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-111

MW: Calculated MW: 50 kDa; Observed MW: 50 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NUP50 in HeLa, A549, HL-60, U251 lysates using NUP50 antibody.



Immunocytochemistry analysis of NUP50 (green) in Jurkat using NUP50 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
