

OGA RABBIT MAB

Cat.#: N262514

Product Name: Anti-OGA Rabbit Monoclonal Antibody

Synonyms: HEXC; MEA5; Mgea5; NCOAT; OGA

UNIPROT ID: O60502

Background: Isoform 1: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins. Isoform 3: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins.

Immunogen: A synthetic peptide of human MGEA5

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-1K5

MW: Calculated MW: 103 kDa; Observed MW: 130 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MGEA5 in HeLa lysates using OGA antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
