

## ORP150 RABBIT MAB

**Cat.#:** N262644

**Product Name:** Anti-ORP150 Rabbit Monoclonal Antibody

**Synonyms:** Grp170; HSP12A; Hyou1; Orp150

**UNIPROT ID:** Q9Y4L1

**Background:** Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular chaperone and participate in protein folding.

**Immunogen:** A synthetic peptide of human ORP150

**Applications:** WB,IHC-P

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R04-4K7

**MW:** Calculated MW: 111 kDa; Observed MW: 150 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cardiovascular

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ORP150 in 293T, MCF-7 lysates using ORP150 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---