

PCSK9 RABBIT MAB

Cat.#: N262682

Product Name: Anti-PCSK9 Rabbit Monoclonal Antibody

Synonyms: HCHOLA3; LDLCQ1; NARC1; Neural apoptosis regulated convertase 1; PC9; PCSK9; Proprotein convertase 9

UNIPROT ID: Q8NBP7

Background: May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.

Immunogen: A synthetic peptide of human PCSK9

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-4F7

MW: Calculated MW: 74 kDa; Observed MW: 65,80 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PCSK9 in HeLa lysates using PCSK9 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
