

PHOSPHO-CSK (SER364) RABBIT PAB

Cat.#: N225328

Product Name: Anti-Phospho-CSK (Ser364) Rabbit pAb

Synonyms: CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL

UNIPROT ID: P41240

Background: Specifically phosphorylates 'Tyr-504' on LCK, which acts as a negative regulatory site. Can also act on the LYN and FYN kinases.

Immunogen: Synthetic peptide of human CSK

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 51 kDa; Observed MW: 55 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Phosphorylated

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-CSK (Ser364) in HUVEC lysates using Phospho-CSK (Ser364) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
