

PHOSPHO-GARNL1 (THR736) RABBIT PAB

Cat.#: N225506

Product Name: Anti-Phospho-GARNL1 (Thr736) Rabbit pAb

Synonyms: p240; GRIPE; GARNL1; TULIP1; RalGAPalpha1

UNIPROT ID: Q6GYQ0

Background: Catalytic subunit of the heterodimeric RalGAP1 complex which acts as a GTPase activator for the Ras-like small GTPases RALA and RALB.

Immunogen: Synthesized peptide derived from human GARNL1 (Phospho-T736)

Applications: IHC-P, WB

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 230 kDa; Observed MW: 220 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human, Mouse, Rat

Conjugation: Unconjugated

Modification: Phosphorylated

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing