

PKNOX2 RABBIT PAB

Cat.#: N225462

Product Name: Anti-PKNOX2 Rabbit pAb

Synonyms: PKNOX2; PREP2; Homeobox protein PKNOX2; Homeobox protein PREP-2; PBX/knotted homeobox 2

UNIPROT ID: Q96KN3

Background: Homeodomain proteins are sequence-specific transcription factors that share a highly conserved DNA-binding domain and play fundamental roles in cell proliferation, differentiation, and death.

Immunogen: Synthesized peptide derived from the C-terminal region of human PREP-2. at AA range: 310-390

Applications: WB,IHC-P,ELISA

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 52 kDa; Observed MW: 55 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PKNOX2 in HepG2 lysates using PKNOX2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
