

PPCEL RABBIT MAB

Cat.#: N263502

Product Name: Anti-PPCEL Rabbit Monoclonal Antibody

Synonyms: CMS22; Prolylendopeptidase-like

UNIPROT ID: Q4J6C6

Background: Serine peptidase whose precise substrate specificity remains unclear.

Immunogen: Recombinant protein of human PPCEL/PREPL

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-8A2

MW: Calculated MW: 84 kDa; Observed MW: 72 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat,Hamster


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PPCEL/PREPL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PPCEL/PREPL in mouse brain lysates using PPCEL antibody.



Western blot analysis of PPCEL/PREPL in HeLa, CHO-K1, C6 lysates using PPCEL/PREPL antibody.