

PPCS RABBIT MAB

Cat.#: N263504

Product Name: Anti-PPCS Rabbit Monoclonal Antibody

Synonyms: CMD2C

UNIPROT ID: Q9HAB8

Background: Catalyzes the second step in the biosynthesis of coenzyme A from vitamin B5, where cysteine is conjugated to 4'-phosphopantothenate to form 4-phosphopantothenoylcysteine.

Immunogen: A synthetic peptide of human PPCS

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-6I3

MW: Calculated MW: 34 kDa; Observed MW: 34 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PPCS in Hela, Jurkat lysates using PPCS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
