

PREALBUMIN RABBIT MAB

Cat.#: N262755

Product Name: Anti-Prealbumin Rabbit Monoclonal Antibody

Synonyms: Prealbumin; TBPA; Amyloidosis I; Transthyretin

UNIPROT ID: P02766

Background: Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Immunogen: A synthetic peptide of human Prealbumin

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-7D9

MW: Calculated MW: 16 kDa; Observed MW: 16 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Prealbumin in Human tonsil, Human placenta lysates using Prealbumin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
